meRgreeMS framework TMT Run 1 TMT Run 1 TMT Run 2 TMT Run 2 TMT Run 1 TMT Run 2 phosphorylated unmodified lysin acetylated phosphorylated unmodified lysin acetylated .xlsx .xlsx .xlsx .xlsx .xlsx .xlsx Remove PSMs: 1. with identification q-value < 0.01 **Filtering** 2. without reporter ions labels ("No Quan Labels") pdf 3. having PTMs of ptmRS localization score < 90 4. that the reference channel is NA PSMs QC Report Perform for every sample: Normalization 1. log2 transformation pdf and 2. zero center normalization by median Rescale 3. rescale to the common reference channel PSMs QC Report Perform for each report file: PSM to Peptide 1. peptide sequences retrieval from PD sequences Summarization 2. re-annotation of PTMs 3. averaging PSM ratios per peptide Peptides QC . Report Perform for each run: Aggregation of Peptides Aggregation aggregation of peptides of Peptides found in multiple files as for Run 1 for Run 2 the average of ratios Aggregation pdf of Runs Peptides QC Report .csv Merged dataset